

CDH8 Rabbit Polyclonal Antibody

Catalog No: #54593

Package Size: #54593-1 50ul #54593-2 100ul



SAB[®]
Signalway Antibody


Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description	
Product Name	CDH8 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant protein of human CDH8.
Conjugates	Unconjugated
Other Names	—
Accession No.	Swiss Prot:_Gene ID:_
Calculated MW	—
SDS-PAGE MW	88kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB□1:500 - 1:2000IF□1:50 - 1:200


Images



Western blot analysis of extracts of various cell lines, using CDH8 Rabbit pAb.



Immunofluorescence analysis of C6 cells using CDH8 Rabbit pAb.



Immunofluorescence analysis of NIH3T3 cells using CDH8 Rabbit pAb.

Background

—

Note: This product is for in vitro research use only and is not intended for use in humans or animals.