

FAM160B2 Rabbit Polyclonal Antibody

Catalog No: #54598

Package Size: #54598-1 50ul #54598-2 100ul


Description

Product Name	FAM160B2 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of Mouse FAM160B2.
Conjugates	Unconjugated
Other Names	G430067P06Rik, Rai16
Accession No.	Swiss Prot:Q80YR2Gene ID:239170
Calculated MW	—
SDS-PAGE MW	82kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB□1:500 - 1:2000IHC□1:50 - 1:200IF□1:50 - 1:200


Images



Western blot analysis of extracts of various cell lines, using FAM160B2 antibody.



Immunohistochemistry of paraffin-embedded mouse kidney using FAM160B2 Rabbit pAb.



Immunofluorescence analysis of NIH-3T3 cells using FAM160B2 Rabbit pAb.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.