

TRRAP Rabbit Polyclonal Antibody

Catalog No: #54600

Package Size: #54600-1 50ul #54600-2 100ul


Description

Product Name	TRRAP Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Rat
Immunogen Description	Recombinant protein of human TRRAP.
Conjugates	Unconjugated
Other Names	PAF350/400, PAF400, STAF40, TR-AP, Tra1
Accession No.	Swiss Prot:Q9Y4A5Gene ID:8295
Calculated MW	—
SDS-PAGE MW	434kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB□1:2000 - 1:5000IHC□1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using TRRAP antibody.



Immunohistochemistry of paraffin-embedded human esophageal cancer using TRRAP Rabbit pAb.

Background

This gene encodes a large multidomain protein of the phosphoinositide 3-kinase-related kinases (PIKK) family. The encoded protein is a common component of many histone acetyltransferase (HAT) complexes and plays a role in transcription and DNA repair by recruiting HAT complexes to chromatin. Deregulation of this gene may play a role in several types of cancer including glioblastoma multiforme. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Sep 2011]

Note: This product is for in vitro research use only and is not intended for use in humans or animals.