

CHML Rabbit Polyclonal Antibody

Catalog No: #54614

Package Size: #54614-1 50ul #54614-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CHML Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant protein of human CHML.
Conjugates	Unconjugated
Other Names	REP2;CHML
Accession No.	Swiss Prot:P26374Gene ID:1122
Calculated MW	—
SDS-PAGE MW	Refer to figures
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB□1:500 - 1:2000IF□1:50 - 1:200

Images



Immunofluorescence analysis of L929 cells using CHML antibody.

Background

—

Note: This product is for in vitro research use only and is not intended for use in humans or animals.