

GRID2 Rabbit Polyclonal Antibody

Catalog No: #54619

Package Size: #54619-1 50ul #54619-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GRID2 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human GRID2.
Conjugates	Unconjugated
Other Names	GluD2;SCAR18;GRID2
Accession No.	Swiss Prot:O43424Gene ID:2895
Calculated MW	—
SDS-PAGE MW	Refer to figures
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB□1:500 - 1:2000IF□1:50 - 1:200


Images



Immunofluorescence analysis of C6 cells using GRID2 antibody.



Immunofluorescence analysis of L929 cells using GRID2 antibody.



Immunofluorescence analysis of U-2 OS cells using GRID2 antibody.

Background

—

Note: This product is for in vitro research use only and is not intended for use in humans or animals.