

NDUFAB1 Rabbit Polyclonal Antibody

Catalog No: #54621

Package Size: #54621-1 50ul #54621-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NDUFAB1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant protein of human NDUFAB1.
Conjugates	Unconjugated
Other Names	NDUFAB1;ACP;ACP1;FASN2A;SDAP
Accession No.	Swiss Prot:O14561Gene ID:4706
Calculated MW	17Da
SDS-PAGE MW	Refer to figures
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200

Images



Immunofluorescence analysis of L929 cells using NDUFAB1 antibody.

Background

-

Note: This product is for in vitro research use only and is not intended for use in humans or animals.