

POU4F2 Rabbit Polyclonal Antibody

Catalog No: #54624

Package Size: #54624-1 50ul #54624-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	POU4F2 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human POU4F2.
Conjugates	Unconjugated
Other Names	POU4F2;BRN3.2;BRN3B;Brn-3b
Accession No.	Swiss Prot:Q12837Gene ID:5458
Calculated MW	29kDa/43kDa
SDS-PAGE MW	Refer to figures
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200

Images



Immunofluorescence analysis of U-2 OS cells using POU4F2 antibody.

Background

The protein encoded by this gene is a member of the POU-domain transcription factor family and may be involved in maintaining visual system neurons in the retina. The level of the encoded protein is also elevated in a majority of breast cancers, resulting in accelerated tumor growth.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.