

ZNF207 Rabbit Polyclonal Antibody

Catalog No: #54646

Package Size: #54646-1 50ul #54646-2 100ul


Description

Product Name	ZNF207 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human ZNF207.
Conjugates	Unconjugated
Other Names	BuGZ; hBuGZ
Accession No.	Swiss Prot:O43670Gene ID:7756
Calculated MW	51kDa
SDS-PAGE MW	51KDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:2000 - 1:5000IHC 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using ZNF207 antibody.

Background

—

Note: This product is for in vitro research use only and is not intended for use in humans or animals.