

HNRNPA3 Rabbit Polyclonal Antibody

Catalog No: #54652

Package Size: #54652-1 50ul #54652-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	HNRNPA3 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic Peptide of human HNRNPA3
Conjugates	Unconjugated
Other Names	HNRNPA3;2610510D13Rik;D10S102;FBRNP;HNRPA3
Accession No.	Swiss Prot:P51991Gene ID:220988
Calculated MW	37kDa/39kDa
SDS-PAGE MW	40KDa
Formulation	Buffer: PBS with 0.02% sodium azide,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB□1:2000 - 1:5000IHC□1:50 - 1:200


Images




Western blot analysis of extracts of U-251MG cells, using HNRNPA3 antibody.



Immunohistochemistry of paraffin-embedded rat ovary using HNRNPA3 Rabbit pAb.



Immunohistochemistry of paraffin-embedded human lung cancer using HNRNPA3 Rabbit pAb.



Immunohistochemistry of paraffin-embedded mouse brain using HNRNPA3 Rabbit pAb.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.