

Fam50b Rabbit Polyclonal Antibody

Catalog No: #54653

Package Size: #54653-1 50ul #54653-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Fam50b Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of Mouse Fam50b.
Conjugates	Unconjugated
Other Names	X5; X5L; XAP-5-; D0H6S2654; D0H6S2654E; XAP-5-like
Accession No.	Swiss Prot:Q9WTJ8Gene ID:108161
Calculated MW	40kDa
SDS-PAGE MW	39KDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:2000 - 1:5000IHC 1:50 - 1:200


Images




Western blot analysis of extracts of Mouse testis, using Fam50b antibody.




Western blot analysis of extracts of Rat testis, using Fam50b antibody.



Immunohistochemistry of paraffin-embedded human esophageal cancer using Fam50b Rabbit pAb.



Immunohistochemistry of paraffin-embedded mouse lung using Fam50b Rabbit pAb.



Immunohistochemistry of paraffin-embedded rat lung using Fam50b Rabbit pAb.

Background

—

Note: This product is for in vitro research use only and is not intended for use in humans or animals.