

EIF4A1/EIF4A2/EIF4A3 Rabbit Polyclonal Antibody

Catalog No: #54656

Package Size: #54656-1 50ul #54656-2 100ul


Description

Product Name	EIF4A1/EIF4A2/EIF4A3 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human EIF4A1/EIF4A2/EIF4A3.
Conjugates	Unconjugated
Other Names	—
Accession No.	Swiss Prot:P60842/Q14240/P38919Gene ID:1973/1974/9775
Calculated MW	—
SDS-PAGE MW	46KDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200


Images




Western blot analysis of extracts of various cell lines, using EIF4A1\EIF4A2\EIF4A3 antibody.



Immunofluorescence analysis of C6 cells using
EIF4A1\EIF4A2\EIF4A3 Rabbit pAb.



Immunofluorescence analysis of NIH-3T3 cells using
EIF4A1\EIF4A2\EIF4A3 Rabbit pAb.



Immunofluorescence analysis of U-2 OS cells using
EIF4A1\EIF4A2\EIF4A3 Rabbit pAb.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.