

FADD Rabbit Polyclonal Antibody

Catalog No: #54657

Package Size: #54657-1 50ul #54657-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FADD Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human FADD.
Conjugates	Unconjugated
Other Names	Mort1/F; Mort1/FADD
Accession No.	Swiss Prot:Q61160Gene ID:14082
Calculated MW	—
SDS-PAGE MW	28KDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB□1:500 - 1:2000IF□1:50 - 1:200

Images



Western blot analysis of extracts of Mouse heart, using FADD antibody.



Western blot analysis of extracts of Jurkat cells, using FADD antibody.

Background

—

Note: This product is for in vitro research use only and is not intended for use in humans or animals.