

BOP1 Rabbit Polyclonal Antibody

Catalog No: #54659

Package Size: #54659-1 50ul #54659-2 100ul


Description

Product Name	BOP1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human BOP1 (NP_056016.1).
Conjugates	Unconjugated
Other Names	—
Accession No.	Swiss Prot:Q14137Gene ID:23246
Calculated MW	84kDa
SDS-PAGE MW	110KDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB□1:2000 - 1:5000IHC□1:50 - 1:200


Images



Western blot analysis of extracts of various cell lines, using BOP1 antibody.



Immunohistochemistry of paraffin-embedded rat lung using BOP1 Rabbit pAb.



Immunohistochemistry of paraffin-embedded mouse brain using BOP1 Rabbit pAb.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.