

RPS6KA1/RPS6KA2/RPS6KA3 Rabbit Polyclonal Antibody

Catalog No: #54663

Package Size: #54663-1 50ul #54663-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	RPS6KA1/RPS6KA2/RPS6KA3 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human RPS6KA1/RPS6KA2/RPS6KA3.
Conjugates	Unconjugated
Other Names	—
Accession No.	Swiss Prot:Q15418/Q15349/P51812Gene ID:6195/ 6196/ 6197
Calculated MW	—
SDS-PAGE MW	83KDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB□1:500 - 1:2000IHC□1:50 - 1:200IF□1:50 - 1:200


Images




Western blot analysis of extracts of various cell lines, using RPS6KA1\RPS6KA2\RPS6KA3 antibody.




Immunohistochemistry of paraffin-embedded human liver using RPS6KA1\RPS6KA2\RPS6KA3 Rabbit pAb.



Immunohistochemistry of paraffin-embedded mouse kidney using RPS6KA1\RPS6KA2\RPS6KA3 Rabbit pAb.



Immunohistochemistry of paraffin-embedded rat stomach using RPS6KA1\RPS6KA2\RPS6KA3 Rabbit pAb.



Immunofluorescence analysis of NIH-3T3 cells using RPS6KA1\RPS6KA2\RPS6KA3 Rabbit pAb.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.