

RAB11A/RAB11B Rabbit Polyclonal Antibody

Catalog No: #54664

Package Size: #54664-1 50ul #54664-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RAB11A/RAB11B Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human RAB11A/RAB11B.
Conjugates	Unconjugated
Other Names	—
Accession No.	Swiss Prot:P62491/Q15907Gene ID:8766/ 9230
Calculated MW	—
SDS-PAGE MW	25KDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200IF 1:50 - 1:200


Images




Western blot analysis of extracts of various cell lines, using RAB11A/RAB11B antibody.




Western blot analysis of extracts of various cell lines, using RAB11A/RAB11B antibody.




Immunohistochemistry of paraffin-embedded human brain using RAB11A/RAB11B Rabbit pAb.



Immunohistochemistry of paraffin-embedded human liver using RAB11A/RAB11B Rabbit pAb.



Immunohistochemistry of paraffin-embedded mouse heart using RAB11A/RAB11B Rabbit pAb.



Immunohistochemistry of paraffin-embedded rat brain using RAB11A/RAB11B Rabbit pAb.

—

Note: This product is for in vitro research use only and is not intended for use in humans or animals.