

SMYD4 Rabbit Polyclonal Antibody

Catalog No: #54679

Package Size: #54679-1 50ul #54679-2 100ul


Description

Product Name	SMYD4 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human SMYD4 (NP_443160.2).
Conjugates	Unconjugated
Other Names	SMYD4;ZMYND21
Accession No.	Swiss Prot:Q8IYR2Gene ID:114826
Calculated MW	89kDa
SDS-PAGE MW	90kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using SMYD4 antibody.



Immunofluorescence analysis of MCF-7 cells using SMYD4 antibody.

Background

-

Note: This product is for in vitro research use only and is not intended for use in humans or animals.