

DVL3 Rabbit Polyclonal Antibody

Catalog No: #54799

Package Size: #54799-1 50ul #54799-2 100ul


Description

Product Name	DVL3 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human DVL3
Conjugates	Unconjugated
Other Names	DVL3;DRS3
Accession No.	Uniprot:Q92997GenelD:1857
Calculated MW	76kDa/78kDa
SDS-PAGE MW	93KDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:1000IHC 1:50 - 1:100IF 1:20 - 1:50


Images




Western blot analysis of extracts of various cell lines, using DVL3 antibody.




Immunohistochemistry of paraffin-embedded mouse brain using DVL3 Rabbit pAb.



Immunohistochemistry of paraffin-embedded mouse kidney using DVL3 Rabbit pAb.



Immunohistochemistry of paraffin-embedded rat ovary using DVL3 Rabbit pAb.



Immunohistochemistry of paraffin-embedded rat brain using DVL3 Rabbit pAb.

Background

This gene is a member of a multi-gene family which shares strong similarity with the *Drosophila* dishevelled gene, *dsh*. The *Drosophila* dishevelled gene encodes a cytoplasmic phosphoprotein that regulates cell proliferation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.