

RAB5C Rabbit Polyclonal Antibody

Catalog No: #55189

Package Size: #55189-1 50ul #55189-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RAB5C Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human RAB5C (NP_958842.1).
Conjugates	Unconjugated
Other Names	RAB5C;L1880;RAB5CL;RAB5L;RABL
Accession No.	Uniprot:P51148GenelD:5878
Calculated MW	23kDa/27kDa
SDS-PAGE MW	23kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using RAB5C antibody.

Background

Members of the Rab protein family are small GTPases of the Ras superfamily that are thought to ensure fidelity in the process of docking and/or fusion of vesicles with their correct acceptor compartment (Han et al., 1996 [PubMed 8646882]).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.