

LGR6 Rabbit Polyclonal Antibody

Catalog No: #55200

Package Size: #55200-1 50ul #55200-2 100ul


Description

Product Name	LGR6 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human LGR6 (NP_067649.2).
Conjugates	Unconjugated
Other Names	LGR6;GPCR;VTS20631
Accession No.	Uniprot:Q9HBX8GeneID:59352
Calculated MW	89kDa/99kDa/104kDa
SDS-PAGE MW	104kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using LGR6 antibody.

Background

This gene encodes a member of the leucine-rich repeat-containing subgroup of the G protein-coupled 7-transmembrane protein superfamily. The encoded protein is a glycoprotein hormone receptor with a large N-terminal extracellular domain that contains leucine-rich repeats important for the formation of a horseshoe-shaped interaction motif for ligand binding. Alternative splicing of this gene results in multiple transcript variants.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.