

ATG10 Rabbit Polyclonal Antibody

Catalog No: #55203

Package Size: #55203-1 50ul #55203-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ATG10 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human ATG10 (NP_113670.1).
Conjugates	Unconjugated
Other Names	ATG10;APG10;APG10L;pp12616
Accession No.	Swiss Prot:Q9H0Y0GeneID:83734
Calculated MW	14kDa/25kDa
SDS-PAGE MW	25kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB□1:500 - 1:2000IHC□1:50 - 1:200

Images



Western blot analysis of extracts of mouse small intestine, using ATG10 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.