

PDCD2L Rabbit Polyclonal Antibody

Catalog No: #55205

Package Size: #55205-1 50ul #55205-2 100ul


Description

Product Name	PDCD2L Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human PDCD2L (NP_115722.1).
Conjugates	Unconjugated
Other Names	PDCD2L
Accession No.	Swiss Prot:Q9BRP1GeneID:84306
Calculated MW	39kDa
SDS-PAGE MW	45kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200


Images



Western blot analysis of extracts of various cell lines, using PDCD2L at 1:500 dilution.



Immunohistochemistry of paraffin-embedded rat heart using PDCD2L at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human liver using PDCD2L at dilution of 1:100 (40x lens).

Background

—

Note: This product is for in vitro research use only and is not intended for use in humans or animals.