

RAB5B Rabbit Polyclonal Antibody

Catalog No: #55221

Package Size: #55221-1 50ul #55221-2 100ul


Description

Product Name	RAB5B Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human RAB5B (NP_002859.1).
Conjugates	Unconjugated
Other Names	RAB5B
Accession No.	Swiss Prot:P61020GenelD:5869
Calculated MW	19kDa/23kDa
SDS-PAGE MW	24kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200IF 1:50 - 1:200


Images



Western blot analysis of extracts of various cell lines, using RAB5B at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit (RM00021). Exposure time: 90s.



Immunofluorescence analysis of L929 cells using RAB5B at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of human colon carcinoma using RAB5B at dilution of 1:100. Blue: DAPI for nuclear staining.

Background

—

Note: This product is for in vitro research use only and is not intended for use in humans or animals.