

RPS6KL1 Rabbit Polyclonal Antibody

Catalog No: #55231

Package Size: #55231-1 50ul #55231-2 100ul


Description

Product Name	RPS6KL1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human RPS6KL1 (NP_113652.2).
Conjugates	Unconjugated
Other Names	RPS6KL1;RSKL2
Accession No.	Uniprot:Q9Y6S9GeneID:83694
Calculated MW	28kDa/59kDa/60kDa
SDS-PAGE MW	60KDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200

Images



Western blot analysis of extracts of Rat brain cells, using RPS6KL1 antibody.



Immunofluorescence analysis of A549 cells using RPS6KL1 antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.