

RAC3 Rabbit Polyclonal Antibody

Catalog No: #55236

Package Size: #55236-1 50ul #55236-2 100ul


Description

Product Name	RAC3 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human RAC3 (NP_005043.1).
Conjugates	Unconjugated
Other Names	RAC3
Accession No.	Uniprot:P60763GenelD:5881
Calculated MW	21kDa
SDS-PAGE MW	21kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using RAC3 antibody.

Background

The protein encoded by this gene is a GTPase which belongs to the RAS superfamily of small GTP-binding proteins. Members of this superfamily appear to regulate a diverse array of cellular events, including the control of cell growth, cytoskeletal reorganization, and the activation of protein kinases. Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.