

SESN2 Rabbit Polyclonal Antibody

Catalog No: #55243

Package Size: #55243-1 50ul #55243-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SESN2 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human SESN2 (NP_113647.1).
Conjugates	Unconjugated
Other Names	SESN2;HI95;SES2;SEST2;sestrin-2
Accession No.	Uniprot:P58004GenelD:83667
Calculated MW	54kDa
SDS-PAGE MW	57kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using SESN2 antibody.

Background

This gene encodes a member of the sestrin family of PA26-related proteins. The encoded protein may function in the regulation of cell growth and survival. This protein may be involved in cellular response to different stress conditions.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.