

COX1 Rabbit Polyclonal Antibody

Catalog No: #55247

Package Size: #55247-1 50ul #55247-2 100ul


Description

Product Name	COX1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	A synthetic peptide of mouse COX1 (NP_904330.1).
Conjugates	Unconjugated
Other Names	MT-CO1;COI;MTCO1;COX1
Accession No.	Uniprot:/GeneID:17708
Calculated MW	56kDa
SDS-PAGE MW	44kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of mouse skeletal muscle, using MT-CO1 antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.