

CBR3 Rabbit Polyclonal Antibody

Catalog No: #55249

Package Size: #55249-1 50ul #55249-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CBR3 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human CBR3 (NP_001227.1).
Conjugates	Unconjugated
Other Names	CBR3;HEL-S-25;SDR21C2;hCBR3
Accession No.	Uniprot:O75828GeneID:874
Calculated MW	30kDa
SDS-PAGE MW	31kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using CBR3 antibody.

Background

Carbonyl reductase 3 catalyzes the reduction of a large number of biologically and pharmacologically active carbonyl compounds to their corresponding alcohols. The enzyme is classified as a monomeric NADPH-dependent oxidoreductase. CBR3 contains three exons spanning 11.2 kilobases and is closely linked to another carbonyl reductase gene - CBR1.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.