

INHBC Rabbit Polyclonal Antibody

Catalog No: #55254

Package Size: #55254-1 50ul #55254-2 100ul


Description

Product Name	INHBC Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human INHBC (NP_005529.1).
Conjugates	Unconjugated
Other Names	INHBC;IHBC
Accession No.	Uniprot:P55103GenelD:3626
Calculated MW	38kDa
SDS-PAGE MW	38kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using INHBC antibody.

Background

This gene encodes a member of the TGF-beta (transforming growth factor-beta) superfamily of proteins. The encoded preproprotein is proteolytically processed to generate a subunit of homodimeric and heterodimeric activin complexes. The heterodimeric complex may function in the inhibition of activin A signaling. Transgenic mice overexpressing this gene exhibit defects in testis, liver and prostate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.