

SH3GLB2 Rabbit Polyclonal Antibody

Catalog No: #55259

Package Size: #55259-1 50ul #55259-2 100ul


Description

Product Name	SH3GLB2 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human SH3GLB2 (NP_064530.1).
Conjugates	Unconjugated
Other Names	SH3GLB2;PP6569;PP9455;RRIG1
Accession No.	Uniprot:Q9NR46GeneID:56904
Calculated MW	43kDa/44kDa
SDS-PAGE MW	44kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using SH3GLB2 antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.