

CENPQ Rabbit Polyclonal Antibody

Catalog No: #55264

Package Size: #55264-1 50ul #55264-2 100ul


Description

Product Name	CENPQ Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human CENPQ (NP_060602.2).
Conjugates	Unconjugated
Other Names	CENPQ;C6orf139;CENP-Q
Accession No.	Uniprot:Q7L2Z9GeneID:55166
Calculated MW	30kDa
SDS-PAGE MW	31kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of Jurkat cells, using CENPQ antibody.

Background

CENPQ is a subunit of a CENPH (MIM 605607)-CENPI (MIM 300065)-associated centromeric complex that targets CENPA (MIM 117139) to centromeres and is required for proper kinetochore function and mitotic progression (Okada et al., 2006 [PubMed 16622420]).[supplied by OMIM, Mar 2008]

Note: This product is for in vitro research use only and is not intended for use in humans or animals.