

IGSF3 Rabbit Polyclonal Antibody

Catalog No: #55285

Package Size: #55285-1 50ul #55285-2 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	IGSF3 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human IGSF3 (NP_001007238.1).
Conjugates	Unconjugated
Other Names	IGSF3;EWI-3;LCDD;V8
Accession No.	Uniprot:O75054GeneID:3321
Calculated MW	135kDa/137kDa
SDS-PAGE MW	125kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using IGSF3 Antibody.

Background

The protein encoded by this gene is an immunoglobulin-like membrane protein containing several V-type Ig-like domains. A mutation in this gene has been associated with bilateral nasolacrimal duct obstruction (LCDD).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.