

CLIP1 Rabbit Polyclonal Antibody

Catalog No: #55298

Package Size: #55298-1 50ul #55298-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CLIP1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human CLIP1 (NP_002947.1).
Conjugates	Unconjugated
Other Names	CLIP1;CLIP;CLIP-170;CLIP170;CYLN1;RSN
Accession No.	Uniprot:P30622GenelD:6249
Calculated MW	156kDa/160kDa/162kDa
SDS-PAGE MW	162kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using CLIP1 antibody.

Background

The protein encoded by this gene links endocytic vesicles to microtubules. This gene is highly expressed in Reed-Sternberg cells of Hodgkin disease. Several transcript variants encoding different isoforms have been found for this gene.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.