

PIP5K1B Rabbit Polyclonal Antibody

Catalog No: #55313

Package Size: #55313-1 50ul #55313-2 100ul


Description

Product Name	PIP5K1B Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human PIP5K1B (NP_003549.1).
Conjugates	Unconjugated
Other Names	PIP5K1B;MSS4;STM7
Accession No.	Swiss Prot:O14986GeneID:8395
Calculated MW	56kDa/61kDa/62kDa
SDS-PAGE MW	71kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using PIP5K1B at 1:1000 dilution.



Immunofluorescence analysis of U2OS cells using PIP5K1B at dilution of 1:100. Blue: DAPI for nuclear staining.

Background

—

Note: This product is for in vitro research use only and is not intended for use in humans or animals.