

COG4 Rabbit Polyclonal Antibody

Catalog No: #55329

Package Size: #55329-1 50ul #55329-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	COG4 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human COG4 (NP_056201.2).
Conjugates	Unconjugated
Other Names	COG4;CDG2J;COD1
Accession No.	Uniprot:Q9H9E3GeneID:25839
Calculated MW	38kDa/81kDa/89kDa
SDS-PAGE MW	89kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of mouse testis, using COG4 antibody.

Background

The protein encoded by this gene is a component of an oligomeric protein complex involved in the structure and function of the Golgi apparatus. Defects in this gene may be a cause of congenital disorder of glycosylation type IIj. Two transcript variants encoding different isoforms have been found for this gene.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.