

RAB17 Rabbit Polyclonal Antibody

Catalog No: #55339

Package Size: #55339-1 50ul #55339-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RAB17 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human RAB17 (NP_071894.1).
Conjugates	Unconjugated
Other Names	RAB17
Accession No.	Uniprot:Q9H0T7GeneID:64284
Calculated MW	9kDa/23kDa
SDS-PAGE MW	23kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using RAB17 antibody.

Background

The Rab subfamily of small GTPases plays an important role in the regulation of membrane trafficking. RAB17 is an epithelial cell-specific GTPase (Lutcke et al., 1993 [PubMed 8486736]).[supplied by OMIM, Oct 2009]

Note: This product is for in vitro research use only and is not intended for use in humans or animals.