

HSPA2 Rabbit Polyclonal Antibody

Catalog No: #55366

Package Size: #55366-1 50ul #55366-2 100ul


Description

Product Name	HSPA2 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human HSPA2 (NP_068814.2).
Conjugates	Unconjugated
Other Names	HSPA2;HSP70-2;HSP70-3
Accession No.	Swiss Prot:P54652GenelD:3306
Calculated MW	70kDa
SDS-PAGE MW	75kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200


Images



Western blot analysis of extracts of various cell lines, using HSPA2 at 1:1000 dilution.



Immunofluorescence analysis of L929 cells using HSPA2 at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of U2OS cells using HSPA2 at dilution of 1:100. Blue: DAPI for nuclear staining.

Background

—

Note: This product is for in vitro research use only and is not intended for use in humans or animals.