

IFNA6 Rabbit Polyclonal Antibody

Catalog No: #55367

Package Size: #55367-1 50ul #55367-2 100ul


Description

Product Name	IFNA6 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human IFNA6 (NP_066282.1).
Conjugates	Unconjugated
Other Names	IFNA6;IFN-alphaK
Accession No.	Swiss Prot:P05013GenelD:3443
Calculated MW	22kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB□1:500 - 1:2000IF□1:50 - 1:200

Images



Immunofluorescence analysis of U2OS cells using IFNA6 at dilution of 1:100. Blue: DAPI for nuclear staining.

Background

—

Note: This product is for in vitro research use only and is not intended for use in humans or animals.