

CAPRIN1 Rabbit Polyclonal Antibody

Catalog No: #55370

Package Size: #55370-1 50ul #55370-2 100ul


Description

Product Name	CAPRIN1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human CAPRIN1
Conjugates	Unconjugated
Other Names	CAPRIN1;GPIAP1;GPIP137;GRIP137;M11S1;RNG105;p137GPI;caprin-1
Accession No.	Swiss Prot:Q14444GeneID:4076
Calculated MW	76kDa/78kDa
SDS-PAGE MW	120kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200IF 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using CAPRIN1 at 1:500 dilution.



Immunohistochemistry of paraffin-embedded mouse brain using CAPRIN1 (40x lens).

Background

-

Note: This product is for in vitro research use only and is not intended for use in humans or animals.