

UCHL5 Rabbit Polyclonal Antibody

Catalog No: #55392

Package Size: #55392-1 50ul #55392-2 100ul


Description

Product Name	UCHL5 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human UCHL5 (NP_001186190.1).
Conjugates	Unconjugated
Other Names	UCHL5;CGI-70;INO80R;UCH-L5;UCH37
Accession No.	Swiss Prot:Q9Y5K5GeneID:51377
Calculated MW	36kDa/37kDa
SDS-PAGE MW	38kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:100

Images



Western blot analysis of extracts of various cell lines, using UCHL5 at 1:1000 dilution.



Immunofluorescence analysis of A549 cells using UCHL5 .
Blue: DAPI for nuclear staining.

Background

—

Note: This product is for in vitro research use only and is not intended for use in humans or animals.