

MT-ND5 Rabbit Polyclonal Antibody

Catalog No: #55436

Package Size: #55436-1 50ul #55436-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MT-ND5 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of mouse MT-ND5 (NP_904338.1).
Conjugates	Unconjugated
Other Names	MT-ND5;MTND5;ND5
Accession No.	Swiss Prot:P03915GenelD:17721
Calculated MW	68kDa
SDS-PAGE MW	67kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB□1:500 - 1:2000IHC□1:50 - 1:200


Images




Western blot analysis of extracts of various cell lines, using MT-ND5 at 1:1000 dilution.




Immunohistochemistry of paraffin-embedded rat brain using MT-ND5 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded rat kidney using MT-ND5 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human colon carcinoma using MT-ND5 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human prostate cancer using MT-ND5 at dilution of 1:100 (40x lens).

Background

—

Note: This product is for in vitro research use only and is not intended for use in humans or animals.