

USF2 Rabbit Polyclonal Antibody

Catalog No: #55441

Package Size: #55441-1 50ul #55441-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	USF2 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human USF2 (NP_003358.1).
Conjugates	Unconjugated
Other Names	USF2;FIP;bHLHb12
Accession No.	Swiss Prot:Q15853GeneID:7392
Calculated MW	23kDa/30kDa/36kDa
SDS-PAGE MW	40kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000 IF □ 1:50 - 1:100

Images



Western blot analysis of extracts of Rat liver, using USF2 at 1:1000 dilution.

Background

This gene encodes a member of the basic helix-loop-helix leucine zipper family of transcription factors. The encoded protein can activate transcription through pyrimidine-rich initiator (Inr) elements and E-box motifs and is involved in regulating multiple cellular processes.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.