

SLC25A11 Rabbit Polyclonal Antibody

Catalog No: #55444

Package Size: #55444-1 50ul #55444-2 100ul


Description

Product Name	SLC25A11 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human SLC25A11 (NP_003553.2).
Conjugates	Unconjugated
Other Names	SLC25A11;OGC;SLC20A4
Accession No.	Swiss Prot:Q02978GeneID:8402
Calculated MW	28kDa/34kDa
SDS-PAGE MW	34kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:100IF 1:50 - 1:100


Images




Western blot analysis of extracts of various cell lines, using SLC25A11 at 1:1000 dilution.




Immunohistochemistry of paraffin-embedded rat spleen using SLC25A11 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human colon carcinoma using SLC25A11 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human gastric cancer using SLC25A11 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse spinal cord using SLC25A11 at dilution of 1:100 (40x lens).

Background

—

Note: This product is for in vitro research use only and is not intended for use in humans or animals.