

THOC1 Rabbit Polyclonal Antibody

Catalog No: #55449

Package Size: #55449-1 50ul #55449-2 100ul


Description

Product Name	THOC1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human THOC1 (NP_005122.2).
Conjugates	Unconjugated
Other Names	THOC1;HPR1;P84;P84N5
Accession No.	Swiss Prot:Q96FV9GeneID:9984
Calculated MW	43kDa/75kDa
SDS-PAGE MW	92kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:100IF 1:50 - 1:100


Images




Western blot analysis of extracts of various cell lines, using THOC1 at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit (RM00021). Exposure time: 30s.




Immunofluorescence analysis of U2OS cells using THOC1 at dilution of 1:100. Blue: DAPI for nuclear staining.




Immunohistochemistry of paraffin-embedded mouse kidney using THOC1 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded rat kidney using THOC1 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human lung cancer using THOC1 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human stomach using THOC1 at dilution of 1:100 (40x lens).

—

Note: This product is for in vitro research use only and is not intended for use in humans or animals.