

METTL7A Rabbit Polyclonal Antibody

Catalog No: #55451

Package Size: #55451-1 50ul #55451-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	METTL7A Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human METTL7A (NP_054752.3).
Conjugates	Unconjugated
Other Names	METTL7A;AAM-B
Accession No.	Swiss Prot:Q9H8H3GeneID:25840
Calculated MW	28kDa
SDS-PAGE MW	28kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB□1:500 - 1:2000IHC□1:50 - 1:200IF□1:50 - 1:200


Images




Immunohistochemistry of paraffin-embedded rat liver using METTL7A at dilution of 1:200 (40x lens).



Immunohistochemistry of paraffin-embedded human prostate using METTL7A at dilution of 1:200 (40x lens).



Immunohistochemistry of paraffin-embedded human stomach using METTL7A at dilution of 1:200 (40x lens).



Western blot analysis of extracts of 22Rv1 cells, using METTL7A at 1:1000 dilution.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.