

## COPS7A Rabbit Polyclonal Antibody

Catalog No: #55454

Package Size: #55454-1 50ul #55454-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |                                                           |
|-----------------------|-----------------------------------------------------------|
| Product Name          | COPS7A Rabbit Polyclonal Antibody                         |
| Host Species          | Rabbit                                                    |
| Clonality             | Polyclonal                                                |
| Isotype               | IgG                                                       |
| Purification          | Affinity purification                                     |
| Applications          | WB                                                        |
| Species Reactivity    | Human,Mouse,Rat                                           |
| Immunogen Description | Recombinant fusion protein of human COPS7A (NP_057403.1). |
| Conjugates            | Unconjugated                                              |
| Other Names           | COPS7A;CSN7;CSN7A;SGN7a                                   |
| Accession No.         | Uniprot:Q9UBW8GenelD:50813                                |
| Calculated MW         | 30kDa                                                     |
| SDS-PAGE MW           | 30kDa                                                     |
| Formulation           | PBS with 0.02% sodium azide,50% glycerol,pH7.3.           |
| Storage               | Store at -20°C. Avoid freeze / thaw cycles.               |

## Application Details

WB 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using COPS7A antibody.

## Background

This gene encodes a component of the COP9 signalosome, an evolutionarily conserved multi-subunit protease that regulates the activity of the ubiquitin conjugation pathway. Alternatively spliced transcript variants that encode the same protein have been described.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.