

HAUS4 Rabbit Polyclonal Antibody

Catalog No: #55457

Package Size: #55457-1 50ul #55457-2 100ul


Description

Product Name	HAUS4 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human HAUS4 (NP_060285.2).
Conjugates	Unconjugated
Other Names	HAUS4;C14orf94
Accession No.	Uniprot:Q9H6D7GeneID:54930
Calculated MW	27kDa/28kDa/37kDa/42kDa
SDS-PAGE MW	36kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using HAUS4 antibody.

Background

This gene encodes a subunit of the centrosome complex termed the human augmin complex. The encoded protein localizes to the spindle microtubules and may play a role in mitotic spindle assembly and maintenance of centrosome integrity during cell division. Alternate splicing results in multiple transcript variants. A pseudogene of this gene is found on chromosome 1.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.