

CEP55 Rabbit Polyclonal Antibody

Catalog No: #55458

Package Size: #55458-1 50ul #55458-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CEP55 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human CEP55 (NP_060601.3).
Conjugates	Unconjugated
Other Names	CEP55;C10orf3;CT111;URCC6
Accession No.	Swiss Prot:Q53EZ4GeneID:55165
Calculated MW	46kDa/54kDa
SDS-PAGE MW	60kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB□1:500 - 1:2000IF□1:50 - 1:200

Images



Western blot analysis of extracts of THP-1 cells, using CEP55 at 1:1000 dilution.



Immunofluorescence analysis of L929 cells using CEP55 at dilution of 1:100. Blue: DAPI for nuclear staining.

Background

–

Note: This product is for in vitro research use only and is not intended for use in humans or animals.