

TMBIM1 Rabbit Polyclonal Antibody

Catalog No: #55462

Package Size: #55462-1 50ul #55462-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TMBIM1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse
Immunogen Description	A synthetic peptide of human TMBIM1 (NP_071435.2).
Conjugates	Unconjugated
Other Names	TMBIM1;LFG3;MST100;MSTP100;PP1201;RECS1
Accession No.	Swiss Prot:Q969X1Gene ID:64114
Calculated MW	34kDa
SDS-PAGE MW	35kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB□1:500 - 1:2000IF□1:50 - 1:200

Images



Western blot analysis of extracts of U-251MG cells, using RECS1 antibody.



Immunofluorescence analysis of NIH-3T3 cells using TMBIM1 Rabbit pAb.

Background

—

Note: This product is for in vitro research use only and is not intended for use in humans or animals.