

# GLYATL1 Rabbit Polyclonal Antibody

Catalog No: #55467

Package Size: #55467-1 50ul #55467-2 100ul


## Description

|                       |                                                               |
|-----------------------|---------------------------------------------------------------|
| Product Name          | GLYATL1 Rabbit Polyclonal Antibody                            |
| Host Species          | Rabbit                                                        |
| Clonality             | Polyclonal                                                    |
| Isotype               | IgG                                                           |
| Purification          | Affinity purification                                         |
| Applications          | WB,IF                                                         |
| Species Reactivity    | Human,Mouse                                                   |
| Immunogen Description | Recombinant fusion protein of human GLYATL1 (NP_001207423.1). |
| Conjugates            | Unconjugated                                                  |
| Other Names           | GLYATL1;GATF-C;GNAT                                           |
| Accession No.         | Swiss Prot:Q969I3Gene ID:92292                                |
| Calculated MW         | 35kDa/38kDa                                                   |
| SDS-PAGE MW           | 35kDa                                                         |
| Formulation           | Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.       |
| Storage               | Store at -20°C. Avoid freeze / thaw cycles.                   |


## Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200

## Images



Western blot analysis of extracts of mouse liver, using GLYATL1 antibody.



Immunofluorescence analysis of U2OS cells using GLYATL1 antibody.

## Background

—

Note: This product is for in vitro research use only and is not intended for use in humans or animals.