

SHC4 Rabbit Polyclonal Antibody

Catalog No: #55471

Package Size: #55471-1 50ul #55471-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SHC4 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human SHC4 (NP_976224.3).
Conjugates	Unconjugated
Other Names	SHC4;RaLP;SHCD
Accession No.	Swiss Prot:Q6S5L8Gene ID:399694
Calculated MW	42kDa/68kDa
SDS-PAGE MW	Refer to figures
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB□1:500 - 1:2000IF□1:50 - 1:200

Images



Immunofluorescence analysis of U2OS cells using SHC4 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.